



EST1A rabbit pAb

Catalog No	YP-Ab-09781
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SMG6 C17orf31 EST1A KIAA0732
Protein Name	Telomerase-binding protein EST1A (EC 3.1.-.-) (EST1-like protein A) (Ever shorter telomeres 1A) (Smg-6 homolog) (Telomerase subunit EST1A) (hSmg5/7a)
Immunogen	Synthesized peptide derived from human EST1A AA range: 688-738
Specificity	This antibody detects endogenous levels of human EST1A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus, nucleolus . Chromosome, telomere . Cytoplasm, cytosol . Particularly enriched in the nucleolus. .
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
