



TPK1 rabbit pAb

Catalog No	YP-Ab-10002
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TPK1
Protein Name	Thiamin pyrophosphokinase 1 (hTPK1) (EC 2.7.6.2) (Placental protein 20) (PP20) (Thiamine pyrophosphokinase 1)
Immunogen	Synthesized peptide derived from human TPK1 AA range: 171-221
Specificity	This antibody detects endogenous levels of human TPK1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	Detected in heart, kidney, testis, small intestine and peripheral blood leukocytes, and at very low levels in a variety of tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

