



# VOPP1 rabbit pAb

<b>Catalog No</b>	YP-Ab-10017
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse;Rat
<b>Applications</b>	WB
<b>Gene Name</b>	VOPP1 ECOP GASP
<b>Protein Name</b>	Vesicular, overexpressed in cancer, prosurvival protein 1 (EGFR-coamplified and overexpressed protein) (ECop) (Glioblastoma-amplified secreted protein) (Putative NF-kappa-B-activating protein 055N)
<b>Immunogen</b>	Synthesized peptide derived from human VOPP1 AA range: 120-170
<b>Specificity</b>	This antibody detects endogenous levels of human VOPP1
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	Cytoplasmic vesicle membrane ; Single-pass type I membrane protein . When overexpressed, localizes in the nucleus and perinuclear regions.
<b>Tissue Specificity</b>	Widely expressed with highest levels in thymus and ovary.
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

---

---

---

---

---

---

---

---