



MSL2 rabbit pAb

Catalog No	YP-Ab-10036
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MSL2 KIAA1585 MSL2L1 RNF184
Protein Name	Male-specific lethal 2 homolog (MSL-2) (Male-specific lethal 2-like 1) (MSL2-like 1) (Male-specific lethal-2 homolog 1) (RING finger protein 184)
Immunogen	Synthesized peptide derived from human MSL2 AA range: 9-59
Specificity	This antibody detects endogenous levels of human MSL2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

