



ANKH1 rabbit pAb

Catalog No	YP-Ab-10046
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ANKHD1 KIAA1085 MASK VBARP PP2500
Protein Name	Ankyrin repeat and KH domain-containing protein 1 (HIV-1 Vpr-binding ankyrin repeat protein) (Multiple ankyrin repeats single KH domain) (hMASK)
Immunogen	Synthesized peptide derived from human ANKH1 AA range: 435-485
Specificity	This antibody detects endogenous levels of human ANKH1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm .
Tissue Specificity	Ubiquitous with high expression in cervix, spleen and brain. Expressed in hematopoietic cells with increased expression in leukemia cells. Isoform 2 is highly expressed in spleen with almost no expression in muscle and brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

