



PKHG3 rabbit pAb

Catalog No	YP-Ab-10094
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PLEKHG3 KIAA0599
Protein Name	Pleckstrin homology domain-containing family G member 3 (PH domain-containing family G member 3)
Immunogen	Synthesized peptide derived from human PKHG3 AA range: 557-607
Specificity	This antibody detects endogenous levels of human PKHG3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm, cytoskeleton . Colocalizes with actin at the leading edge of polarized cells. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

