



MSTP9 rabbit pAb

Catalog No	YP-Ab-10127
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	MST1P9 D1F15S1A MSTP9
Protein Name	Putative macrophage-stimulating protein MSTP9 (Brain rescue factor 1) (BRF-1) (Hepatocyte growth factor-like protein homolog)
Immunogen	Synthesized peptide derived from human MSTP9 AA range: 154-204
Specificity	This antibody detects endogenous levels of human MSTP9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

