



TMEM9 rabbit pAb

Catalog No	YP-Ab-10140
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TMEM9 DERP4 HSPC186 PSEC0012 UNQ631/PRO1248
Protein Name	Transmembrane protein 9 (Dermal papilla-derived protein 4)
Immunogen	Synthesized peptide derived from human TMEM9 AA range: 29-79
Specificity	This antibody detects endogenous levels of human TMEM9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Lysosome membrane ; Single-pass type I membrane protein . Late endosome membrane ; Single-pass type I membrane protein . Endosome, multivesicular body membrane ; Single-pass type I membrane protein .
Tissue Specificity	Highly expressed in adrenal gland, thyroid gland, testis, ovary and prostate (PubMed:12359240). Moderate expression in trachea, spinal cord, stomach, colon, small intestine and spleen (PubMed:12359240, PubMed:30374053). Low expression in bone marrow, lymph node, thymus and peripheral blood lymphocytes (PubMed:12359240). Expression is detected in hematopoietic cell lines including those of myeloid, erythroid, B- and T-cell origin (PubMed:12359240).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images