



ZNHI2 rabbit pAb

Catalog No	YP-Ab-10142
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ZNHIT2 C11orf5
Protein Name	Zinc finger HIT domain-containing protein 2 (Protein FON)
Immunogen	Synthesized peptide derived from human ZNHI2 AA range: 26-76
Specificity	This antibody detects endogenous levels of human ZNHI2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	Low expression in most tissues; highly expressed in testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

