



SPOT1 rabbit pAb

Catalog No	YP-Ab-10147
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SPANXA2-OT1 CXorf18
Protein Name	Putative uncharacterized protein SPANXA2-OT1
Immunogen	Synthesized peptide derived from human SPOT1 AA range: 11-61
Specificity	This antibody detects endogenous levels of human SPOT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

