



CS057 rabbit pAb

Catalog No	YP-Ab-10156
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	C19orf57
Protein Name	Uncharacterized protein C19orf57 (Pre-T/NK cell-associated protein 3B3)
Immunogen	Synthesized peptide derived from human CS057 AA range: 246-296
Specificity	This antibody detects endogenous levels of human CS057
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Chromosome . During meiosis, recruited to chromosomes and localizes on recombination sites in a double-strand break-dependent manner. First appears on the chromosome axis at leptotene. Along with the progression of meiotic recombination, released from the axis to form bridge-like structures linking homolog axes before they are synapsed. Finally, located between synapsed homolog axes and on the synaptonemal complex (SC). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images