



U1SBP rabbit pAb

Catalog No	YP-Ab-10171
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SNRNP35 HM1 U1SNRNPBP
Protein Name	U11/U12 small nuclear ribonucleoprotein 35 kDa protein (U11/U12 snRNP 35 kDa protein) (U11/U12-35K) (Protein HM-1) (U1 snRNP-binding protein homolog)
Immunogen	Synthesized peptide derived from human U1SBP AA range: 168-218
Specificity	This antibody detects endogenous levels of human U1SBP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in heart, liver, skeletal muscle and pancreas.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

