



TM158 rabbit pAb

Catalog No	YP-Ab-10189
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMEM158 HBBP RIS1
Protein Name	Transmembrane protein 158 (40 kDa BINP-binding protein) (p40BBP) (Ras-induced senescence protein 1)
Immunogen	Synthesized peptide derived from human TM158 AA range: 97-147
Specificity	This antibody detects endogenous levels of human TM158
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
