



SDIM1 rabbit pAb

Catalog No	YP-Ab-10208
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SDIM1
Protein Name	Stress-responsive DNAJB4-interacting membrane protein 1
Immunogen	Synthesized peptide derived from human SDIM1 AA range: 83-133
Specificity	This antibody detects endogenous levels of human SDIM1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in brain with higher detection in neurons than astrocytes. Decreased expression in Alzheimer brains. Detected at protein level in brain and cervix.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

