



## ZKSC3 rabbit pAb

<b>Catalog No</b>	YP-Ab-10242
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	ZKSCAN3 ZFP47 ZNF306 ZNF309 ZSCAN13
<b>Protein Name</b>	Zinc finger protein with KRAB and SCAN domains 3 (Zinc finger and SCAN domain-containing protein 13) (Zinc finger protein 306) (Zinc finger protein 309) (Zinc finger protein 47 homolog) (Zf47) (Zfp-47)
<b>Immunogen</b>	Synthesized peptide derived from human ZKSC3 AA range: 457-507
<b>Specificity</b>	This antibody detects endogenous levels of human ZKSC3
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	Nucleus. Cytoplasm. Mainly localizes in the nucleus. Under starvation conditions translocates to the cytoplasm, allowing expression of target genes involved in autophagy and lysosome biogenesis/function.
<b>Tissue Specificity</b>	
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images