



BEST3 rabbit pAb

Catalog No	YP-Ab-10253
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	BEST3 VMD2L3
Protein Name	Bestrophin-3 (Vitelliform macular dystrophy 2-like protein 3)
Immunogen	Synthesized peptide derived from human BEST3 AA range: 79-129
Specificity	This antibody detects endogenous levels of human BEST3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Present in skeletal muscle and weakly in brain, spinal cord, bone marrow and retina.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

