




# Tubulin alpha rabbit pAb(Caenorhabditis elegans)

<b>Catalog No</b>	YP-Ab-10299
<b>Isotype</b>	IgG
<b>Reactivity</b>	C.elegans
<b>Applications</b>	WB
<b>Gene Name</b>	
<b>Protein Name</b>	Tubulin alpha
<b>Immunogen</b>	Purified recombinant Full lenth Tubulin alpha (Caenorhabditis elegans)
<b>Specificity</b>	This antibody detects endogenous levels of Tubulin alpha (Caenorhabditis elegans)
<b>Formulation</b>	Liquid in PBS containing 50% glycerol and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:5000-20000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	Alpha Tubulin
<b>Observed Band</b>	
<b>Cell Pathway</b>	
<b>Tissue Specificity</b>	
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images



Western Blot analysis of C.elegans cells using Tubulin alpha rabbit pAb(Caenorhabditis elegans). Secondary antibody(catalog#:RS0002) was diluted at 1:20000