



Thymidine kinase Polyclona Antibody

Catalog No	YP-Ab-10879
Isotype	IgG
Reactivity	Human virus
Applications	ELISA
Gene Name	TK
Protein Name	Thymidine kinase
Immunogen	Synthesized peptide derived from human Thymidine kinase
Specificity	This antibody detects endogenous levels of human Thymidine kinase
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	ELISA(peptide) 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

