



ANR22 rabbit pAb

Catalog No	YP-Ab-11114
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ANKRD22
Protein Name	ANR22
Immunogen	Synthesized peptide derived from human ANR22 AA range: 65-115
Specificity	This antibody detects endogenous levels of ANR22 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 4 ANK repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western blot analysis of lysates from HEK293 cells,
primary antibody was diluted at 1:1000, 4° over night

