



## RN112 rabbit pAb

<b>Catalog No</b>	YP-Ab-11164
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human; Mouse;Rat
<b>Applications</b>	WB;IHC
<b>Gene Name</b>	RNF112 BFP ZNF179
<b>Protein Name</b>	RN112
<b>Immunogen</b>	Synthesized peptide derived from human RN112 AA range: 377-427
<b>Specificity</b>	This antibody detects endogenous levels of RN112 at Human/Mouse/Rat
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:500-2000;IHC-p 1:50-300
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	Membrane ; Multi-pass membrane protein . Membrane ; Peripheral membrane protein . Cytoplasm . Nucleus . Nucleus, nuclear body . Nucleus, nucleoplasm . Endosome . Cytoplasmic vesicle, secretory vesicle, synaptic vesicle . Cell junction, synapse, postsynaptic density . Perikaryon . Cell projection, neuron projection . Predominantly in the nucleus, but some amounts were also found in the cytoplasm. Oxidative stress stimulates its shuttling from the cytoplasm into the nucleus. Recruited to nuclear bodies via its interaction with ZBTB16. Localizes to the cell soma and neuritis and only slightly to the nucleus in the neurons of most brain areas. .
<b>Tissue Specificity</b>	Predominantly expressed in brain (PubMed:10574464). Decreased expression in glioma brain tumors as compared to normal brains (at protein level) (PubMed:28684796).
<b>Function</b>	similarity:Contains 1 RING-type zinc finger.,tissue specificity:Predominantly expressed in brain.,
<b>Background</b>	This gene encodes a member of the RING finger protein family of transcription factors. The protein is primarily expressed in brain. The gene is located within the Smith-Magenis syndrome region on chromosome 17. [provided by RefSeq, Jul 2008],




matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).