



FAT3 rabbit pAb

Catalog No	YP-Ab-11741
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	IHC;IF
Gene Name	FAT3 CDHF15 KIAA1989
Protein Name	FAT3
Immunogen	Synthesized peptide derived from human FAT3 AA range: 896-946
Specificity	This antibody detects endogenous levels of FAT3 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1: 50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Expressed in ES cells, primitive neuroectoderm, fetal brain, infant brain, adult neural tissues and prostate.
Function	function:May play a role in the interactions between neurites derived from specific subsets of neurons during development.,similarity:Contains 1 laminin G-like domain.,similarity:Contains 33 cadherin domains.,similarity:Contains 4 EGF-like domains.,tissue specificity:Expressed in ES cells, primitive neuroectoderm, fetal brain, infant brain, adult neural tissues and prostate.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images