




## F179B rabbit pAb

<b>Catalog No</b>	YP-Ab-12086
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human; Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	FAM179B KIAA0423
<b>Protein Name</b>	F179B
<b>Immunogen</b>	Synthesized peptide derived from human F179B AA range: 1664-1714
<b>Specificity</b>	This antibody detects endogenous levels of F179B at Human/Mouse
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1: 500-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	Cell projection, cilium . Cytoplasm, cytoskeleton . Cytoplasm, cytoskeleton, cilium axoneme . Detected along the length of primary cilia and at the basal body. Colocalization with the cytoplasmic microtubule cytoskeleton upon heterologous expression is most likely an artifact. .
<b>Tissue Specificity</b>	
<b>Function</b>	similarity:Belongs to the FAM179 family.,similarity:Contains 6 HEAT repeats.,
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images



Western blot analysis of lysates from MCF-7 cells,  
primary antibody was diluted at 1:1000, 4° over night