




DANCR rabbit pAb

Catalog No	YP-Ab-12253
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	DANCR ANCR KIAA0114 SNHG13
Protein Name	DANCR
Immunogen	Synthesized peptide derived from human DANCR AA range: 52-102
Specificity	This antibody detects endogenous levels of DANCR at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	Expressed in keratinocytes.
Function	caution:Product of a dubious CDS prediction.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4° over night