




BMP4 rabbit pAb

Catalog No	YP-Ab-12465
Isotype	IgG
Reactivity	Mouse;Rat
Applications	WB; ELISA
Gene Name	BMP4
Protein Name	BMP4
Immunogen	Synthesized peptide derived from human BMP4
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat BMP4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Effects of bone morphogenetic protein 4 on TGF- β 1-induced cell proliferation, apoptosis, activation and differentiation in mouse lung fibroblasts via ERK/p38 MAPK signaling pathway
PeerJ
Wang Wenju Lu WB Mouse Jian