



Ankyrin B rabbit pAb

Catalog No	YP-Ab-12540
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ANK2
Protein Name	Ankyrin B
Immunogen	Synthesized peptide derived from human Ankyrin B
Specificity	This antibody detects endogenous levels of Ankyrin B at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.24% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ankyrin-2 (ANK-2) (Ankyrin-B) (Brain ankyrin) (Non-erythroid ankyrin)
Observed Band	
Cell Pathway	Cytoplasm, cytoskeleton . Membrane . Cytoplasm, myofibril, sarcomere, M line . Apical cell membrane . Cell membrane . Cell junction, synapse, postsynaptic cell membrane . Early endosome . Recycling endosome . Lysosome . Mitochondrion . Cytoplasm, myofibril, sarcomere, Z line . Cell membrane, sarcolemma, T-tubule . Expressed at the apical membrane of airway lung epithelial cells (By similarity). Localized to the plasma membrane of the inner segments of photoreceptors in retina. Colocalizes with SPTBN1 in a distinct intracellular compartment of neonatal cardiomyocytes (PubMed:19007774). In skeletal muscle, localizes to neuromuscular junctions (By similarity). Localizes with puncta at mitochondria ends. Colocalizes and cotransports on motile vesicles with RABGAP1L (By similarity). .
Tissue Specificity	Present in plasma membrane of neurons as well as glial cells throughout the brain. Expressed in fetal brain and in temporal cortex of adult brain. Also expressed in the inner segments of rod photoreceptors in retina.
Function	
Background	



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images