



CHCH1 rabbit pAb

Catalog No	YP-Ab-17269
Isotype	IgG
Reactivity	Human, Mouse
Applications	IHC, WB
Gene Name	CHCHD1 C10orf34
Protein Name	Coiled-coil-helix-coiled-coil-helix domain-containing protein 1 (Nuclear protein C2360)
Immunogen	Synthesized peptide derived from human N-terminal CHCH1
Specificity	This antibody detects endogenous levels of CHCH1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit, polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Coiled-coil-helix-coiled-coil-helix domain-containing protein 1 (Nuclear protein C2360)
Observed Band	
Cell Pathway	Mitochondrion . Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

