



## Anti-Covid-19 N protein mAb IgM (S06FC)

<b>Catalog No</b>	YP-Ab-17296
<b>Isotype</b>	IgG
<b>Reactivity</b>	Covid-19
<b>Applications</b>	ELISA ECL Immunogold
<b>Gene Name</b>	SARS-CoV-2 N, 2019-nCoV N, Covid-19 N
<b>Protein Name</b>	SARS-CoV-2 Nucleoprotein , 2019-nCoV Nucleoprotein, Covid-19 Nucleoprotein
<b>Immunogen</b>	The recombinant SARS-CoV-2 (2019-nCoV) N Protein
<b>Specificity</b>	This antibody detects SARS-CoV-2 (2019-nCoV) Nucleoprotein
<b>Formulation</b>	PBS
<b>Source</b>	Humanized mAb IgM
<b>Purification</b>	The antibody was affinity-purified from 293T expression by affinity-chromatography using specific immunogen.
<b>Dilution</b>	ELISA 1:5000-1:100000, WB 1:1000-1:10000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	
<b>Tissue Specificity</b>	
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

---

---

---

---

---

---

---

---