



peptide YY rabbit pAb

货号	YP-Ab-17760
同位型	IgG
靶点	PYY
基因名称	PYY
种属	Human;Mouse;Rat
应用	WB
蛋白名称	Peptide YY (PYY) (PYY-I) (Peptide tyrosine tyrosine) [Cleaved into: Peptide YY(3-36) (PYY-II)]
免疫原	Synthesized peptide derived from human peptide YY
特异性	This antibody detects endogenous levels of peptide YY at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	11kD
功能	This gut peptide inhibits exocrine pancreatic secretion, has a vasoconstrictory action and inhibits jejunal and colonic mobility.
细胞定位	Secreted.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images