

## Goat Anti-Rabbit IgG(H+L)(Cyanine3 conjugated)

Catalog No. YPE800085

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

### Product Information

**Immunogen:** Rabbit IgG

**Purification Method:** Affinity purification

**Conjugation:** CY3

**Concentration:** 0.5mg/mL(specific antibody purification $\geq$ 99%),CY3:Ab=3-9:1

**Applications:** IF

**Dilution:** Immunofluorescence: 1:50-100

Optimal working dilutions must be determined by end user. Dilute before when it's needed, and the residual should be discarded.

**Amax:** 554nm

**E<sub>max</sub>:** 568-574nm

**Buffer:** 0.01M PBS with 50% glycerol and 0.05% Sodium Azide PH7.6

**Storage:** Store at -20°C. Avoid repeated freeze / thaw cycles.Protected from prolonged exposure to light.