

Goat Anti-Rabbit IgG(H+L)(FITC conjugated)

Catalog No. YPE800086-F

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Immunogen: Rabbit IgG

Purification Method: Affinity purification

Conjugation: FITC

Concentration: 0.5mg/mL(specific antibody purification \geq 99%),FITC:Ab=4-6:1

Applications: IF

Dilution: Immunofluorescence: 1:20-100

Optimal working dilutions must be determined by end user. Dilute before when it's needed, and the residual should be discarded.

Amax: 495nm

E_{max}: 525nm

Buffer: 0.01M PBS with 50% glycerol and 0.05% Sodium Azide PH7.6

Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.Protected from prolonged exposure to light.