



优品生物

UpingBio

网址: www.upingBio.com

邮箱: youpin@youpinbio.com

地址: 深圳市南山区粤海街道高新区社区高新南九道 45 号

技术电话: 400-999-8863

山羊抗兔免疫球蛋白 G (IgG) H+L (AbFluor 647) 抗体说明书

首先感谢您选择深圳市优品生物科技有限公司，我司为您提供优质的免疫学试剂，在您使用本产品前，请认真阅读说明书，不清楚之处，请联系技术人员解答，这对您顺利完成实验很有帮助。

名称：Goat Anti-Rabbit IgG (H&L) - Alexa Fluor 647

目录号：YPE800092

亚型：IgG

靶点：Rabbit IgG (H&L)

组成：1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

来源：Goat

稀释：IHC (1:200 - 1:1000), IF (1:200 - 1:1000), FCM (1:100 - 1:1000), Elisa (Use at an assay dependent concentration)

纯化工艺：The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

背景：Immunoway secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also available and are recommended for use with immunoglobulin-rich samples.。

应用范围：Elisa;IF;FCM。

种属：Rabbit

产品纯度：1mg/mL。

储存方法：Stable for one year at -15° C to -25° C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap.

Aliquot to avoid repeated freezi.

注意事项: Avoid repeated freezing and melting。

有效期限: Please use it within the validity period of the pipe body。

使用建议: This product can be used for immunological response related experiments, for more information, please consult technicians。