



**优品生物**  
**UpingBio**

网址: [www.upingBio.com](http://www.upingBio.com)

邮箱: [youpin@youpinbio.com](mailto:youpin@youpinbio.com)

地址: 深圳市南山区粤海街道高新区社区高新南九道 45 号

技术电话: 400-999-8863

## 山羊抗人免疫球蛋白G(IgG)抗体 (Fc-AbFluor 488) 说明书

首先感谢您选择深圳市优品生物科技有限公司，我司为您提供优质的免疫学试剂，在您使用本产品前，请认真阅读说明书，不清楚之处，请联系技术人员解答，这对您顺利完成实验很有帮助。

**名称：**山羊抗人免疫球蛋白G(IgG)抗体 (Fc-AbFluor 488)

**目录号：**YPE800097

**亚型：**IgG

**靶点：**Human IgG Fc

**组成：**1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

**来源：**Goat

**稀释：**IHC (1:200 - 1:1000), IF (1:200 - 1:1000), FCM (1:100 - 1:1000), Elisa (Use at an assay dependent concentration)

**纯化工艺：**The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

**其他名称：**Goat Anti-Human IgG Fc - Alexa Fluor 488

**背景：**Immunoway secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also available and are recommended for use with immunoglobulin-rich samples.

**应用范围：**Elisa;IF;FCM

**种属：**Human

**产品纯度：**1mg/mL。

**储存方法:** Stable for one year at -15°C to -25°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezi

**注意事项:** Avoid repeated freezing and melting。

**有效期限:** Please use it within the validity period of the pipe body。

**使用建议:** This product can be used for immunological response related experiments, for more information, please consult technicians。