



**优品生物**  
**UpingBio**

网址: [www.upingBio.com](http://www.upingBio.com)

邮箱: [youpin@youpinbio.com](mailto:youpin@youpinbio.com)

地址: 深圳市南山区粤海街道高新区社区高新南九道 45 号

技术电话: 400-999-8863

## 兔抗羊免疫球蛋白 G (IgG) 抗体 (生物素标记) 说明书

首先感谢您选择深圳市优品生物科技有限公司, 我司为您提供优质的免疫学试剂, 在您使用本产品前, 请认真阅读说明书, 不清楚之处, 请联系技术人员解答, 这对您顺利完成实验很有帮助。

**名称:** 兔抗羊免疫球蛋白 G (IgG) 抗体 (生物素标记)。

**目录号:** YPE800006-S

**亚型:** IgG

**应用:** WB;IHC;IF;ELISA

**种属:** Goat

**靶点:** Goat IgG(h+l)

**组成:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

**来源:** Polyclonal, Rabbit,IgG

**稀释:** WB 1:1000-10000, ELISA 1:10000-20000, IHC 1:50-100 IF 1:50-1000

**纯化工艺:** The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.

**浓度:** 2mg/mL

**储存:** -15° C to -25° C/1 year(Do not lower than -25° C)

**其他名称:** BIOTIN\*Polyclonal Rabbit Anti-Goat IgG(h+l)

**背景:** Immunoway secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also available and are recommended for use with immunoglobulin-rich samples.