



Cleaved-Caspase-4 p20 (Q81) Polyclonal Antibody

Catalog No YP-Ab-00017

Isotype IgG

Reactivity Human

Applications WB;IHC;IF;ELISA

Gene Name

Protein Name

Immunogen

Specificity

Formulation

Source

Purification

Dilution

Concentration 1 mg/ml

Purity ≥90%

Storage Stability -20°C/1 year

Synonyms

Observed Band

Cell Pathway

Tissue Specificity

Function

Background

matters needing attention Avoid repeated freezing and thawing!

Usage suggestions This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

[Products Images](#)