



TBX21 Polyclonal Antibody

Catalog No	YP-Ab-07767
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	TBX21 TBET TBLYM
Protein Name	T-box transcription factor TBX21 (T-box protein 21) (T-cell-specific T-box transcription factor T-bet) (Transcription factor TBLYM)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TBX21 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	58kD
Cell Pathway	Nucleus .
Tissue Specificity	T-cell specific.
Function	disease:Genetic variations in TBX21 are associated with susceptibility to aspirin-induced asthma [MIM:208550]. Asthma is defined as a chronic inflammatory lung disease that is characterized by airway hyperreactivity, eosinophil inflammation, and mucus hypersecretion resulting in intermittent airway obstruction.,disease:Genetic variations in TBX21 are associated with susceptibility to asthma and nasal polyps [MIM:208550]. Nasal polyps is a chronic inflammatory disease of the paranasal sinus mucosa, leading to the protrusion of edematous polyps into the nasal cavities.,function:Transcription factor that controls the expression of the TH1 cytokine, interferon-gamma. Initiates TH1 lineage development from naive TH precursor cells both by activating TH1 genetic programs and by repressing the opposing TH2 programs.,similarity:Contains 1 T-box DNA-binding domain.,tissue specificity:T-cell speci
Background	T-box 21(TBX21) Homo sapiens This gene is a member of a phylogenetically conserved family of genes that share a common DNA-binding domain, the T-box. T-box genes encode transcription factors involved in the regulation of



developmental processes. This gene is the human ortholog of mouse Tbx21/Tbet gene. Studies in mouse show that Tbx21 protein is a Th1 cell-specific transcription factor that controls the expression of the hallmark Th1 cytokine, interferon-gamma (IFNG). Expression of the human ortholog also correlates with IFNG expression in Th1 and natural killer cells, suggesting a role for this gene in initiating Th1 lineage development from naive Th precursor cells. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images