



YB039 rabbit pAb

| | |
|----------------------------------|--|
| Catalog No | YP-Ab-11936 |
| Isotype | IgG |
| Reactivity | Human; Mouse |
| Applications | WB |
| Gene Name | LINC00116 NCRNA00116 |
| Protein Name | YB039 |
| Immunogen | Synthesized peptide derived from human YB039 |
| Specificity | This antibody detects endogenous levels of YB039 at Human/Mouse |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Source | Polyclonal, Rabbit,IgG |
| Purification | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Dilution | WB 1: 500-2000 |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| Synonyms | |
| Observed Band | |
| Cell Pathway | Mitochondrion inner membrane ; Single-pass membrane protein . Preferentially binds to cardiolipin relative to other common cell membrane lipids. . |
| Tissue Specificity | |
| Function | |
| Background | |
| matters needing attention | Avoid repeated freezing and thawing! |
| Usage suggestions | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



Products Images

