



NMBR Polyclonal Antibody

Catalog No	YP-Ab-12629
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;IHC-p;IF(paraffin section)
Gene Name	NMBR
Protein Name	Neuromedin-B receptor
Immunogen	Synthetic Peptide of NMBR
Specificity	The antibody detects endogenous NMBR protein.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1:1000-2000 IHC: 1:200-500
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NMBR; Neuromedin-B receptor; NMB-R; Epididymis tissue protein Li 185a; Neuromedin-B-preferring bombesin receptor; YM3394
Observed Band	50
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in epididymis (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
