



Pdcd-1 rabbit pAb

Catalog No	YP-Ab-15715
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA;IHC-p
Gene Name	PDCD1 PD1
Protein Name	Pdcd-1
Immunogen	Synthesized peptide derived from human Pdcd-1
Specificity	This antibody detects endogenous levels of Human Pdcd-1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Programmed cell death protein 1 (Protein PD-1;hPD-1;CD antigen CD279)
Observed Band	
Cell Pathway	Cell membrane ; Single-pass type I membrane protein.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
