



Antigen Retrieval buffer(20X)TRIS-EDTA of pH8.0

Catalog No YP-Ab-15728

Isotype IgG

Reactivity

Applications IHC

Gene Name

Protein Name

Immunogen

Specificity

Formulation

Source Rabbit

Purification

Dilution Diluted at 1:19 using ddH2O

Concentration 1 mg/ml

Purity ≥90%

Storage Stability -20°C/1 year

Synonyms

Observed Band

Cell Pathway

Tissue Specificity

Function

Background

matters needing attention Avoid repeated freezing and thawing!

Usage suggestions This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

