



CDC37 (phospho-Ser13) rabbit pAb

Catalog No	YP-Ab-16646
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CDC37 CDC37A
Protein Name	CDC37 (Ser13)
Immunogen	Synthesized phospho peptide around human CDC37 (Ser13)
Specificity	This antibody detects endogenous levels of Human Mouse Rat CDC37 (phospho-Ser13)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Hsp90 co-chaperone Cdc37 (Hsp90 chaperone protein kinase-targeting subunit) (p50Cdc37)
Observed Band	44
Cell Pathway	Cytoplasm .
Tissue Specificity	Lymph,Placenta,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
