



Connexin-26 Polyclonal Antibody

Catalog No	YP-Ab-16873
Isotype	IgG
Reactivity	Human;Mouse
Applications	IHC-p;IF(paraffin section)
Gene Name	GJB2
Protein Name	Gap junction beta-2 protein (Connexin-26) (Cx26)
Immunogen	Synthetic Peptide of Connexin-26
Specificity	Connexin-26 protein(A219) detects endogenous levels of Connexin-26
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Dilution	IHC 1:100-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GJB2; Gap junction beta-2 protein; Connexin-26; Cx26; YM3545
Observed Band	26
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Cell junction, gap junction . Colocalizes with GJB4 at gap junction plaques in the cochlea. .
Tissue Specificity	Blood,Colon,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
