



# Rabbit anti-PCYT2 Polyclonal Antibody

|                                  |  |
|----------------------------------|--|
| <b>Catalog No</b>                | YP-Ab-17123  |
| <b>Isotype</b>                   | IgG  |
| <b>Reactivity</b>                | Human,Mouse,Rat  |
| <b>Applications</b>              | WB 1:500 - 1:2000  |
| <b>Gene Name</b>                 | type (类型) : Polyclonal Antibody  |
| <b>Protein Name</b>              |  |
| <b>Immunogen</b>                 | Recombinant fusion protein of human PCYT2 (NP_002852.1).   |
| <b>Specificity</b>               | matters needing attention (注意事项) : Avoid repeated freezing and thawing!  |
| <b>Formulation</b>               |  |
| <b>Source</b>                    | Host (宿主) : Rabbit   |
| <b>Purification</b>              | Affinity purification  |
| <b>Dilution</b>                  |  |
| <b>Concentration</b>             | 1 mg/ml  |
| <b>Purity</b>                    | ≥90%   |
| <b>Storage Stability</b>         | -20°C/1 year   |
| <b>Synonyms</b>                  |  |
| <b>Observed Band</b>             |  |
| <b>Cell Pathway</b>              |  |
| <b>Tissue Specificity</b>        | Usage suggestions (使用建议) : This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.   |
| <b>Function</b>                  | Disclaimer note (免责声明) : The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors. |
| <b>Background</b>                |  |
| <b>matters needing attention</b> | Avoid repeated freezing and thawing!   |
| <b>Usage suggestions</b>         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.  |



## Products Images