



Lambda Light Chain Rabbit pAb

| | |
|----------------------------------|---|
| Catalog No | YP-Ab-17131 |
| Isotype | IgG |
| Reactivity | Human |
| Applications | WB,IHC-P,ICC/IF,FC |
| Gene Name | IGLV1-51/IGLC1 |
| Protein Name | |
| Immunogen | |
| Specificity | |
| Formulation | |
| Source | Rabbit |
| Purification | |
| Dilution | WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100 |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| Synonyms | |
| Observed Band | |
| Cell Pathway | |
| Tissue Specificity | |
| Function | |
| Background | |
| matters needing attention | Avoid repeated freezing and thawing! |
| Usage suggestions | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



Products Images

