



VEGF Receptor 1 Rabbit pAb

Catalog No	YP-Ab-17140
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB,IHC-P,ICC/IF,IP,FC
Gene Name	FLT1
Protein Name	
Immunogen	A synthesized peptide derived from human VEGF Receptor 1
Specificity	
Formulation	
Source	
Purification	Affinity Chromatography
Dilution	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Modification (修改) Unmodified
Observed Band	
Cell Pathway	
Tissue Specificity	Buffer System (缓冲系统) Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Function	
Background	Gene ID: 2321
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

