



EPHX2 rabbit pAb

Catalog No	YP-Ab-17670
Isotype	IgG
Reactivity	EPHX2
Applications	50UL
Gene Name	WB
Protein Name	
Immunogen	EPHX2
Specificity	Bifunctional epoxide hydrolase 2 [Includes: Cytosolic epoxide hydrolase 2 (CEH) (EC 3.3.2.10) (Epoxide hydratase) (Soluble epoxide hydrolase) (SEH); Lipid-phosphate phosphatase (EC 3.1.3.76)] Human Gene Id: 2053 Human Swiss Prot No: P34913 Mouse Gene Id: 13850 Mouse Swiss Prot No: P34914 Rat Gene Id: 65030 Rat Swiss Prot No: P80299
Formulation	Synthesized peptide derived from human EPHX2
Source	This antibody detects endogenous levels of EPHX2 at Human, Mouse,Rat
Purification	Rabbit,polyclonal
Dilution	Liquid in PBS containing 50% glycerol, and 0.569% sodium azide.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	61kD
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images