



PYGB rabbit pAb

货号	YP-Ab-18301
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	PYGB
基因名称	PYGB
蛋白名称	Glycogen phosphorylase, brain form (EC 2.4.1.1)
免疫原	Synthesized peptide derived from human PYGB
特异性	This antibody detects endogenous levels of PYGB at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	93kD
功能	Glycogen phosphorylase that regulates glycogen mobilization . Phosphorylase is an important allosteric enzyme in carbohydrate metabolism . Enzymes from different sources differ in their regulatory mechanisms and in their natural substrates . However, all known phosphorylases share catalytic and structural properties .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images