



MUC20 rabbit pAb

货号	YP-Ab-18302
同位型	IgG
应用	WB
种属	Human,Mouse
宿主	Rabbit
基因名称	MUC20 KIAA1359 UNQ2782/PRO7170
蛋白名称	Mucin-20 (MUC-20)
免疫原	Synthesized peptide derived from human MUC20
特异性	This antibody detects endogenous levels of MUC20 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	78kD
修饰	Unmodified
克隆性	Polyclonal
功能	May regulate MET signaling cascade. Seems to decrease hepatocyte growth factor (HGF)-induced transient MAPK activation. Blocks GRB2 recruitment to MET thus suppressing the GRB2-RAS pathway. Inhibits HGF-induced proliferation of MMP1 and MMP9 expression.
细胞定位	Secreted . Apical cell membrane . Basolateral cell membrane . Cell projection, microvillus membrane .
组织表达	Highly expressed in kidney, moderately in placenta, lung, prostate, liver, and digestive system. In the kidney, localized in the proximal tubules but not in the glomerulus or distal tubules. Detected in most of the male urogenital tract epithelia, with the exception of epididymis.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images