



SNRPF rabbit pAb

货号	YP-Ab-18305
同位型	IgG
应用	WB
种属	Human,Mouse
基因名称	SNRPF PBSCF
蛋白名称	Small nuclear ribonucleoprotein F (snRNP-F) (Sm protein F) (Sm-F) (SmF)
免疫原	Synthesized peptide derived from human SNRPF
特异性	This antibody detects endogenous levels of SNRPF at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	9kD
功能	Plays a role in pre-mRNA splicing as a core component of the spliceosomal U1, U2, U4 and U5 small nuclear ribonucleoproteins (snRNPs), the building blocks of the spliceosome . Component of both the pre-catalytic spliceosome B complex and activated spliceosome C complexes . Is also a component of the minor U12 spliceosome . As part of the U7 snRNP it is involved in histone 3'-end processing .
细胞定位	Cytoplasm, cytosol . Nucleus . SMN-mediated assembly into core snRNPs occurs in the cytosol before SMN-mediated transport to the nucleus to be included in spliceosomes. .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images