



LY6D rabbit pAb

货号	YP-Ab-18307
同位型	IgG
应用	WB
种属	Human,Mouse
靶点	
基因名称	LY6D E48
蛋白名称	Lymphocyte antigen 6D (Ly-6D) (E48 antigen)
免疫原	Synthesized peptide derived from human LY6D
特异性	This antibody detects endogenous levels of LY6D at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	14kD
功能	May act as a specification marker at earliest stage specification of lymphocytes between B- and T-cell development. Marks the earliest stage of B-cell specification.
细胞定位	Cell membrane; Lipid-anchor, GPI-anchor.
组织表达	Expressed exclusively at the outer cell surface of transitional epithelia and the keratinocyte of stratified squamous epithelia.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

