



CFP rabbit pAb

货号	YP-Ab-18308
同位型	IgG
应用	WB
种属	Human,Mouse
靶点	
基因名称	CFP PFC
蛋白名称	Properdin (Complement factor P)
免疫原	Synthesized peptide derived from human CFP
特异性	This antibody detects endogenous levels of CFP at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	52kD
功能	A positive regulator of the alternate pathway (AP) of complement . It binds to and stabilizes the C3- and C5-convertase enzyme complexes . Inhibits CFI-CFH mediated degradation of Complement C3 beta chain (C3b) .
细胞定位	Secreted .
组织表达	
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

